



PELVIPOWER™-TRAINING EQUIPMENT

MAGNETIC FIELD TRAINING

OPTIMAL SEAT POSITION

The individually adjustable seating position provides a comfortable and effective seating position during magnetic field therapy.

TOUCH-SCREEN

The training programs are intuitively set and activated via the easy-to-use touch-screen menu.

ADJUSTABLE ARMREST

The armrest can be folded up if necessary. This facilitates entry, for example, for wheelchair users.

SAFETY

Optional: Automatic detection of metal implants

ADJUSTABLE MAGNETIC FIELD COIL

The magnetic coils located under the seat generate repetitive and short magnetic field impulses. This stimulation reliably strengthens the pelvic floor and the surrounding muscles.

The spool can be moved 18 cm forward and backward. This achieves optimum training for the patient.

